

	Incoming Inspection Specification COOK® Urological		
Title:	Cone Tip		
<small> WARNING - CONFIDENTIAL PROPRIETARY PROPERTY This document is owned by COOK Medical. It contains confidential proprietary trade secret information and must not be copied. The document and the information it contains can be used only by the recipient for the specific use for which it was requested. All other use is strictly prohibited. This document must be returned to COOK Medical immediately upon request by COOK Medical. By possession of this document, the possessor expressly agrees to comply with these terms. "© COPYRIGHT COOK® Urological 2011" </small>	Document Number: CONE TIP	Original Date: 12 December 1990	
	Version Number: 5	Effective Date: 02May2014	
	CR Number: 224-14-T1	Checked By: TI 01May2014	

Name: Cone Tip

Frequency: Refer to [QSI03_05](#) Level II

Reference drawing: Sabin [5932](#) – 5Fr
[5933](#) – 6Fr
[5934](#) – 7Fr
[5935](#) – 8Fr

PROCEDURE:

Inspect all steps 100%, unless otherwise specified. Reference [QSI03_05](#).

1. Using fixture [Y2969](#)-A1, place fixture beneath movable extension. Lift, then lower movable extension until it rests on top of appropriate pin for Fr size cone tip to be checked.

Note: Assembly may need raised or lowered on shaft in order to perform step. 1.
 Loosen screw to raise or lower.

2. Zero dial.
3. Lift movable extension and remove fixture.
4. Place I.D. of cone tip over appropriate pin. Seat.
5. Lift, then lower movable extension until it rest on end of cone tip.
6. Read measurement and compare to print spec.
7. Measure overall length.
8. Check large I.D..
9. Check that parts are the correct color.