

	Controlled Training Procedure Form COOK® Urological		
Title:	Training Guidelines for Bonding Open-End Flexi-Tip Catheters		
<small> WARNING - CONFIDENTIAL PROPRIETARY PROPERTY This document is owned by COOK Medical. It contains confidential proprietary trade secret information and must not be copied. The document and the information it contains can be used only by the recipient for the specific use for which it was requested. All other use is strictly prohibited. This document must be returned to COOK Medical immediately upon request by COOK Medical. By possession of this document, the possessor expressly agrees to comply with these terms. © COPYRIGHT COOK® Urological 2014 </small>	Document Number: CT08-Form_01-Answers	Original Date: 07Feb2014	
	Version Number: 2	Effective Date: 30Jul2015	
	CR Number: CIN-15-3436	Checked By: SD 10Jul2015	

SECTION 1 – WRITTEN EXAM

Trainee: _____

Printed Name	Signature	Date
--------------	-----------	------

Directions: Please answer the following questions. Feel free to use reference specifications to help you answer them.

Trainee must correctly answer the following questions:

1. What Equipment ID numbers can be used to bond 5FR Open-End Flexi-tip Catheters? 1125, 1155, 1172

2. What is used to hold the ID of the tip and tubing to the correct diameter during bonding? Cannula

How is the strength of the tip checked in the catheter department? By flexing the tip

3. True/False The catheter does not need to be butted up against the tip material in the mold.

4. True/~~False~~ The temperature can be changed only within the range specified on the specification.

MI 021305 DWG C0434

Graded by: _____ Overall Result (Pass/Fail)

Printed	Signature	Date
---------	-----------	------
