

	<h1>Controlled Training Procedure Form</h1> <h2>COOK® Urological</h2>	
Title:	Training Guidelines for Bonding Flexi-Tip Catheters on Equip ID 1180	
<p>WARNING - CONFIDENTIAL PROPRIETARY PROPERTY</p> <p>This document is owned by COOK Medical. It contains confidential proprietary trade secret information and must not be copied. The document and the information it contains can be used only by the recipient for the specific use for which it was requested. All other use is strictly prohibited. This document must be returned to COOK Medical immediately upon request by COOK Medical. By possession of this document, the possessor expressly agrees to comply with these terms.</p> <p>© COPYRIGHT COOK® Urological 2013"</p>	<p>Document Number:</p> <p>CT07-Form_01-Answers</p>	<p>Original Date: 20Aug2013</p>
<p>Version Number: 2</p>	<p>Effective Date: 03Aug2015</p>	
<p>CR Number: CIN-15-3436</p>	<p>Checked By: MES 03Aug2015</p>	

SECTION 1 – WRITTEN EXAM

Trainee: _____

Printed Name	Signature	Date
--------------	-----------	------

Directions: Please answer the following questions. Feel free to use reference specifications to help you answer them.

Trainee must correctly answer the following questions:

1. What parameter setting is used for 5fr bonding on Equip ID 1180? 3
2. When changing the heat time what needs to be added to the seconds listed on the specification when entering it into the machine? 0 at then end
3. How are closed end tips loaded into the mold? With tweezers
4. ~~True~~False The catheter does not need to be butted up against the tip material in the mold.
5. ~~True~~False The temperature can be changed only within the range specified on the specification.

MI 02110 SI 428 DWG C4576

Graded by: _____ Overall Result (Pass/Fail)

Printed Signature Date